



CelCradleTM

Cradle for High Cell Density

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

eus@nt-rt.ru || <https://escolifesciences.nt-rt.ru>

The Cradle for Your Cells

CelCradle™ is a cost-effective, single-use, benchtop bioreactor system capable of supporting high density culture of anchorage-dependent or adherent cells. It is a ready-to-use, packed-bed bioreactor system that has linear scalability from laboratory scale to production scale, complete with automated cell harvesting system.

The CelCradle™ operates through the Tide Motion principle wherein cells, attached to BioNOC™ II carriers, are alternately exposed to aeration and nutrition via the decompression and compression of the bellows holding the culture medium. The gentle vertical oscillation of the culture medium creates a dynamic interface between air and culture medium on the surface of the cells, providing the cells with an environment that is of extremely low shear stress, high aeration and nutrition levels, zero foaming, and no O₂ limitation. This efficient nutrient and oxygen transfer allow the CelCradle™ system to produce high cell density and yield.

Applications

The CelCradle™ can be used in many different applications, as exhibited by journal articles (see literature report on page 5). These include the following applications:


- Culture of anchorage-dependent or adherent cells
- Conversion from roller bottles to single-use, closed system
- Overcome limitations of microcarrier-stirred tank bioreactor technology
- Autologous and allogeneic cell therapy
- Human and animal vaccines

5.5 grams of BioNOC™ II (standard)

- Serves as the packed bed where cells adhere to and stabilize within the culture vessel
 - Enhanced biocompatibility
 - Long-lasting hydrophilicity (negatively charged surface for coating)
 - High porosity
 - Low lint content
 - Excellent mechanical strength

Media reservoir

- 500 mL culture medium containing compartment



Dimension: 480 x 300 x 140 mm





Weight: 2 kg

Key Features


- Single-use, pre-sterilized, and ready-to-use bottles
- Compact and small footprint (standard system fits inside a 170 L or 6 ft³ CO₂ incubator)
- Compatible with most serum-free culture medium
- Harvest whole cells or cell components
- Linear scale-up by direct multiplication of bottles or the use of TideXcell®
- Reduced labor, costs, and space requirements
- Large surface area for cell attachment and growth
- Available in batch, fed-batch, and recirculation modes

Established Cell Lines

- | | |
|------------|------------|
| ✓ VERO | ✓ MSC |
| ✓ MDCK | ✓ CHO |
| ✓ MDBK | ✓ XC-18 |
| ✓ Sf9 | ✓ ST |
| ✓ HEK-293 | ✓ CEF |
| ✓ HEK-293T | ✓ GL 37 |
| ✓ BHK 21 | ✓ Marc 145 |



The Tide Motion Principle



Aeration

Bellows are released exposing the carriers to carbon dioxide.
The cells attached to the carriers will take in aeration.



Nutrition

Bellows are pushed up submerging the carriers with culture medium. The cells attached to the carriers will take in nutrition from the culture medium.


Anchorage-dependent cell culture

- One CelCradle™ bottle consists of 5.5 grams of BioNOC™ II providing up to 2,400 cm²/g of surface area for cell attachment and growth. The carriers are proven to have enhanced biocompatibility, long-lasting hydrophilicity, low lint content, and excellent mechanical strength.

Note: surface area may vary depending on the cell line to be cultured.



FDA and Hoescht 33342-stained cells attached to the BioNOC™ II carriers



Autologous and Allogeneic Cell Therapy

Cell Therapy is the use of live whole cells for the treatment of diseases and involves cultivation of cells to a certain density prior to administration. A commonly used cell type in cell therapy is mesenchymal stem cells, which are anchorage-dependent. Currently, cultivation of cells for cell therapy is done using T-flasks; however this method is quite laborious, time and space consuming. T-flasks have limited surface area for growth and thus, require handling of several hundred T-flasks and multiple passaging. The CelCradle™ system's large surface area and compact design help solve these problems, with a single CelCradle™ bottle equalling the productivity of several hundred T-flasks.

Human and Animal Vaccine

Cell culture-based vaccine production is the current trend in vaccine production as it offers several advantages over traditional vaccine production technologies including simpler mass production, rapid manufacturing, independent of SPF animal, controllable quality, and hypoallergenic products.

The CelCradle™ is an ideal system for laboratory-scale, cell-culture based vaccine production as it is capable of supporting high cell density culture, production of high viral titer, and linear scalability to production level. VacciXcell's Tide Motion bioreactor system has been used for the research and development, and production of vaccines for several indications including:

- Influenza
- JEV
- Rabies
- Hog Cholera
- Hepatitis A
- EV71
- Rabies
- IBDV



Apart from these, the CelCradle™ has also been used for other applications including recombinant protein production, pharmacokinetic studies, and cellular component production. Details of these applications can be found from literature support.

Literature Support

The following are some of the available literature support online for the various applications of the CelCradle™ system.


- [1] Asaoka, Y., Tanaka, T., Tsumoto, K., Tomita, M., & Ide, T. (n.d.). Efficient expression of recombinant soluble human FcγRI in mammalian cells and its characterization. *Protein Expression and Purification*, 155-161.
- [2] Brown, A., Singer, D., Mcsharry, J., Barnard, R., Hazuda, D., & Drusano, G. (n.d.). In Vitro Dose Ranging Studies for Serine Protease Inhibitor, MK-4519, Against a Hepatitis C Virus (HCV) Replicon using the BelloCell System. *Antiviral Research*.
- [3] Chen, Y., Wu, J., Wang, K., Chiang, Y., Lai, C., Chung, Y., & Hu, Y. (n.d.). Baculovirus-mediated production of HDV-like particles in BHK cells using a novel oscillating bioreactor. *Journal of Biotechnology*, 135-147.
- [4] Drugmand, J., J.-F., J., Agathos, S., & Schneider, Y. (n.d.). Growth of Mammalian and Lepidopteran Cells on BioNOC® II Disks, a Novel Macroporous Microcarrier. *Cell Technology for Cell Products*, 781-784.
- [5] Hammonds, J., Chen, X., Zhang, X., Lee, F., & Spearman, P. (n.d.). Advances in methods for the production, purification, and characterization of HIV-1 Gag-Env pseudovirion vaccines. *Vaccine*, 8036-8048.
- [6] Ho, L., Greene, C., Schmidt, A., & Huang, L. (n.d.). Cultivation of HEK 293 cell line and production of a member of the superfamily of G-protein coupled receptors for drug discovery applications using a highly efficient novel bioreactor. *Cytotechnology*, 117-123.
- [7] Hu, Y., Lu, J., & Chung, Y. (n.d.). High-density cultivation of insect cells and production of recombinant baculovirus using a novel oscillating bioreactor. *Cytotechnology*, 145-153.
- [8] Huang, K., Lo, W., Chung, Y., Lai, Y., Chen, C., Chou, S., & Hu, Y. (n.d.). Combination of Baculovirus-Mediated Gene Delivery and Packed-Bed Reactor for Scalable Production of Adeno-Associated Virus. *Human Gene Therapy*, 1161-1170.
- [9] Lu, J., Chung, Y., Chan, Z., & Hu, Y. (n.d.). A Novel Oscillating Bioreactor BelloCell: Implications for Insect Cell Culture and Recombinant Protein Production. *Biotechnology Letters Biotechnol Lett*, 1059-1065.
- [10] Mcsharry, J., Singer, D., Kulawy, R., Brown, A., & Drusano, G. (n.d.). Use of the BelloCell System to Determine the Optimal Dose of Ribavirin to Inhibit the Expression of an HCV Replicon in 2209-23 Cells. *Antiviral Research*.
- [11] Toriniwa, H., & Komiya, T. (n.d.). Japanese encephalitis virus production in Vero cells with serum-free medium using a novel oscillating bioreactor. *Biologicals*, 221-226.

Conversion From Roller Bottle to Single-use, Closed System Cell Culture

The roller bottle system plays a major role in cell culture-based vaccine production as it is easy to operate, has a simple scale-up method, observable cell growth and cytopathic effect (CPE), and limited contamination; however, the system has many limitations including intensive labor, large space requirement, high running costs, and low efficiency of culture medium utilization.

In principle, the roller bottle system is very similar to CelCradle™ system. In both systems, cells are alternately exposed to aeration and nutrition. The CelCradle™, however, can overcome the limitations of the roller system while improving output. The advantages of the CelCradle™ system over the roller bottle include:

- One CelCradle™ system is equivalent to eighty roller bottles
- Reduced labor and space requirements
- More controllable
- Reduced contamination risks
- Cell harvest in a closed system
- More efficient use of culture medium
- Higher cell density and viral titer



CelCradle™ System Complete

CelCradle™ System is simple to operate and virtually no learning curve. CelCradle™ 500A bottle is placed in the CelCradle™ unit stage for batch and semi-batch operation where process components are easily traceable.

Filtered Cap

0.22 µm ventilation filter provided in the cap



CelCradle™ Stage

Stainless Steel 316 L Carcass that can hold up to 4 CelCradle™ bottles.

Controller

Magnetized back for convenient positioning on the side or on the front of a CO₂ Incubator.


Place inside a 170 L or 6 ft³ CO₂ incubator



Overcome limitations of microcarrier-stirred tank bioreactor technology

Another technology used for the culture of adherent cells is the use of microcarriers, specially-treated micro-beads where cells attach to and grow on. The microcarriers are suspended in culture medium in stirred tank bioreactors where the medium is continuously agitated and parameters are automatically monitored and controlled. The microcarrier technology, however, has several disadvantages including reduced cell attachment efficiency, shear stress, foaming, and lack of linear scalability, all of which the CelCradle™ system are able to overcome. The CelCradle™ system also has improved features compared to the microcarrier system in terms of the following:

- Viral titer and cell density
- Cell harvest efficiency
- Culture medium usage efficiency
- Harvest purity



BioNOC™ II (packed bed)

Cells remain entrapped in the carriers simplifying media replacement and product harvesting.

Compressible Bellows

Compresses for nutrition and decompresses for aeration.

CelCradle™ Continuous System Complete

Same CelCradle™ unit stage but utilizes CelCradle™ 500AP bottle for recirculation/continuous processes.

Each CelCradle™ System accommodates up to four single-use bottles, making this an ideal screening device to test varying medium formulations or cell lines.

No steam or water lines, autoclave or utilities required, just a power outlet and CO₂ incubator

Fits in a 240 L or 8 ft³ CO₂ incubator

CelCradle™ Stage

Stainless Steel 316 L Carcass that can hold up to 4 CelCradle™ bottles.

Dimension	264 (W) x 359 (L) x 170 (H) mm (10.4 x 14 x 6.7 inches)
Weight	7.0 kg (15.4 lb)
Power	100~230 VOLTS AC, 50/60 Hz (Input); 12 VOLTZ DC (Output)
Up-Down Rate:	0.25 to 2.0 mm/sec. Step of 1.0 sec or 1.0 min.
Delay Time:	0 to 99 min 59 sec Step of 1.0 sec or 1.0 min.
Driver Motor:	DC stepping-motor
Environment	20~42 °C, 0~90% relative humidity (in a CO ₂ incubator)
Mechanical Protection:	Hi-Low Optical Sensor
Transmission:	Gear set (ratio 1:1.2) and belt
Materials:	Aluminum alloy, chromic steel
Recommended CO ₂ incubator dimension	CCL-170L model: 505 x 535 x 633 mm (19.9" x 21.1" x 24.9") and CCL-240L model: 595 x 640 x 633 mm (23.4" x 25.2" x 24.9")

Place inside
a 240 L or 8 ft³
CO₂ incubator



CelFeeder

Enables four-pump operation with individual programming setting.

Media Reservoir

Fresh culture medium bottle in either 1 L or 2 L volume.

Tubing Complete Set (Connection Set)

Connects the media reservoir bottle to the CelCradle™ bottle.

**Disposable Tubing Set and Pump Head**

Replacement tubings and pump head for tubing complete set (connection set).

Culture Medium Sampling Port



Needle closed sampling for cell culture growth monitoring.

Pump Head


Calibration-free autoclavable pump head for up to three times to support consistent feeding rate.




The Main Controller



Message displayed on screen once controller is turned on



Message that will be displayed when stage is rendered on top before opening



Display message after system has warmed up




Dimension	226 (W) x 40 (D) x 137 (L) mm (5.4 x 8.9 x 1.6 inches)
Weight	1.16 kg (2.6 lb)
Power	100~230 VOLTS AC 50/60 Hz (Input) ; 12 VOLTS DC (Output)
Environment	Room temperature (outside CO ₂ incubator)
Materials	Magnetic back plate, to hold the controller to the side of the incubator.

The CelFeeder

CelFeeder is a microprocessor-controlled pumping unit of the CelCradle™ Continuous Complete or other cell culture systems.

The CelFeeder is designed to simplify user's operation such that daily feeding volume and frequency of the feeding are the only two parameters required. It simplifies setup work by avoiding calibration and calculation as usually needed in conventional peristaltic pumps.



- Controls four (4) peristaltic pumps individually
- Controls four (4) recirculation process of CelCradle™ 500P/AP bottles at a time
- Daily feeding rate ranges from 1 – 1999 mL/day
- Feeding frequency of 1-24 times
- Pump calibration made possible to improve accuracy

Dimension	230 (W) x 131 (D) x 61 (H) mm, (5.2 x 9.1 x 2.4 inches)
Weight	0.645 kg (1.42 lb)
Power	85~240 Volts AC, 50/60 Hz (Input); ~180 mA
Environment	Ambient to 45°C, ~95% Relative Humidity



CelFeeder G2™

CelFeeder™ G2 Pump is a HMI-controlled pumping unit for CelCradle™ for different applications. The CelFeeder™ G2 Pump is designed to simplify user's operation such that daily feeding volume and pH/DO monitoring are controlled. It simplifies culture operations, reducing calibration and calculation steps as usually required in conventional peristaltic pumps.

The two front-mounted, fixed-speed peristaltic pumps provide for convenient liquid addition/removal during the culture's recirculation process and automated pH/DO control.

- One CelCradle™ system is equivalent to eighty roller bottles
- Reduced labor and space requirements
- More controllable
- Reduced contamination risks
- More efficient use of culture medium
- Higher cell density and viral titer



Applications:

- pH/DO Control
- Adjusting pH value to optimal conditions
- Recirculation process mode

Features:

1. Sophisticated software
2. Keeps record of the pH/DO data
3. pH control and monitoring
4. Reagent addition for pH control
5. Automated re-circulation of external medium for medium exchange/replenishment during culture
6. Automated pH/DO control of culture environment for optimum cell growth
7. Automatically graphs data in excel-based software
8. 21CFR Part 11 enabled when used with optional Esco Aster Compliant Workbook and interface software package from Ofni Systems (ExcelSafe™).
9. Optional single-use sensor for in line monitoring

SPECIFICATIONS

SPECIFICATIONS				
Dimension		351 (L) x 310 (D) x 400 (H) mm		
Weight		18 Kg		
Power		120 Watts		
Control Software		TIA Portal V16 Siemens		
Working Environment		Ambient temperature		
pH	Length and diameter	S120 mm; 12 mm		
	Range	0-12		
	Precision	± 0.01		
	Probe	Mettler Toledo pH Sensor InPro3253i		
	Control	Liquid addition		
DO	Length and diameter	120 mm; 12 mm		
	Range	0-100%		
	Precision	± 0.01		
	Probe	Mettler Toledo DO Sensor InPro6860i		
	Control	TBD		
Peristaltic Pump	Configurations	Recirculation process or reagent addition		
	Flow Rate	Flow Rates (ml/min) Tubing 1.6 mm wall		
		Tube bore (mm)	Speed (rpm)	
			(ml/rev)	200
		0.5	0.02	4.6
		0.8	0.04	8.6
		1.6	0.14	28
		2.4	0.29	58
		3.2	0.47	95
		4.0	0.67	135
4.8	0.85	170		



Single-Use and Ready-to-Use CelCradle™ Bottles

CelCradle™ bottles are sterilized through gamma irradiation and come pre-packed with 5.5 grams of carriers as standard. Different models of bottles cover a specific application in cell culture.

CelCradle™ 500A



Suitable for cell harvest, transient transfection, slow adhering cells; batch/semi-batch culture

CelCradle™ 500AP



Suitable for cell harvest, transient transfection, slow adhering cells; continuous culture

Application/ Bottle	Item Code	Secreted Protein, Viruses	Cell Harvest (for non-secreted proteins, viruses or cells)	Carrier Harvest (for protein extraction or reuse of carriers)
CelCradle™ 500A	1400003	Applicable, but not optimal	Best application	Best application
CelCradle™ 500AP	1400004	Applicable, but not optimal	Best application	Best application

Cell Monitoring



GlucCell™

A portable, easy-to-use glucose monitoring system for cell culture. Get results in just 15 seconds. Measurement ranges from 30-500 mg/dL.



Crystal Violet Dye

A simple tool for the quantification of cells based on the number of nuclei dyed. The CVD kit contains reagents that disrupt the cells, thereby releasing the nuclei, which are subsequently dyed.

Carrier Sampling

500A/500AP Bottles



Adherent cells attached to the macroporous carriers can be monitored through cell sampling. The bottle is opened inside a biosafety cabinet and using a sterile forceps, carriers can be taken out for cell staining and monitoring.

Inoculation Phase pH Control

OPTION 1

If using a pre-mixed media with fixed amount of sodium bicarbonate (NaHCO_3), readjusting before culture may be difficult. Alternatively, users can opt to:

For Upright Seeding

Seed at a lower CO_2 concentration



For Inverted Seeding

Adjust pH with hydrochloric acid (HCl) to pH 6.8 as a starting point



OPTION 2

Use 500AP bottle and connect one tubing with air filter and the other leaving tightly closed.






attached to a 0.22 μm syringe filter

close cap

OPTION 3

Compress the bellows before closing the bottles with a non-vented cap (white cap). This method lessens excessive air in the CelCradle™ bottle which may affect the pH value of the culture medium. This can be attributed to the decrease of CO₂ concentration in the liquid rendering it basic.



Ordering Information

Item Code	Product Name	Dimension	Weight	Package
2230006	CelCradle™ System Complete	590 x 380 x 420 mm	14 kg	1 x CelCradle™ Stage 1 x GlucCell™ Glucose Monitoring System
2230007	CelCradle™ Continuous System Complete	590 x 380 x 420 mm	15 kg	1 x CelCradle™ Stage 1 x CelFeeder Pump 2 x Tubing Complete Set
2230005	CelCradle™ Stage	360 x 265 x 170 mm	9.88 kg	1x Main Console 1x Control Box 1x 100-240 V Power Adapter 1x Signal Cable 1x Manual CD 2x Forceps 1x Crystal Violet Dye Nucleus Count Kit
1400067	CelFeeder Pump	290 x 210 x 80 mm	1 kg	1 x CelFeeder Pump
1400011	Tubing Complete Set	370 x 150 x 75 mm	0.5 kg	1x Disposable Tubing Accessory 1x Pump Head
1400012	Disposable Tubing Set & Pump Head	370 x 150 x 75 mm	0.5 kg	1x Tubing Set 1x Pump Head
1400013	Disposable Tubing Accessory	370 x 150 x 75 mm	0.5 kg	5x Disposable Tubing Accessory
1400009	GlucCell™ Glucose Monitoring System	96 x 60 x 18.5 mm	0.5 kg	1 x GlucCell™ Glucose meter 2 bottles of Glucose Test Strip (1bt = 25 test strips)
1400010	GlucCell™ Glucose Test Strip	80 x 70 x 35 mm	0.05 kg	2 x Glucose Test Strip bts (2 x 25 pcs)
1400014	Crystal Violet Dye Nucleus Count Kit	50 x 50 x 105 mm	0.5 kg	1 x CVD Bottle (100 ml/bt)
1400015	Filtered Cap	370 x 150 x 75 mm	0.5 kg	Cap for CelCradle™ Bottle
1400016	Non-Vented Cap	370 x 150 x 75 mm	0.5 kg	Cap for CelCradle™ Bottle
1400017	Forceps	250 x 15 x 10 mm	0.1 kg	Used for sampling of BioNOC™ II carriers
1400021	CelCradle™ 500A/500AP Strainer	370 x 150 x 75 mm	0.5 kg	Single-use Strainer (Pack of 10)

CelCradle™ and Esco CO₂ Incubator

The Perfect Combination for High Density Adherent Culture

The CelCradle™ system can be incorporated into an existing CO₂ incubator or can be purchased with an Esco CO₂ incubator. Esco offers a wide range of CO₂ incubators that suit clients' different requirements, provide superior performance and cell protection.

CelCulture® CO₂ Incubator

Sleek, reliable, and intuitive, Esco CelCulture® CO₂ incubators provide all-rounded sample protection that brings your scientific dreams one step closer to reality.

Features:

- HPA-validated 90°C overnight moist heat decontamination cycle
- Rounded corners and seamless design for easy cleaning
- Constructed using electrogalvanized with ISOCIDE™ powder-coating to eliminate 99.9% of surface bacteria within 24 hours of exposure
- With optional copper interiors for added antimicrobial protection
- Available in 50 L, 170 L, and 240 L sizes
- ULPA Filter
 - 99.999% efficient, superior to conventional HEPA filters
 - Filters air continuously
 - Chamber returns to ISO Class 5 cleanliness in 11 minutes upon door closing to prevent contamination
- Direct Heat & Air Jacket
 - Fast and uniform heating
 - Rapid temperature recovery without overshoot
 - Air jacket improves chamber stability
- O₂ Sensor (for suppressed O₂ Model)
 - Long life
 - Stable output signal
 - No influence from CO₂
- CO₂ Sensor
 - CelCradle™ and CO₂ Incubator Combination
 - Single-beam, dual-wavelength IR sensor is drift-free



- Auto-zeroing
- SmartSense™ Microcontroller Interface
 - Intuitive, fully equipped control and monitoring system

CelCulture® CO₂ Incubators IR Sensor Model with Stainless Steel Chamber

230 VAC, 50/60 Hz		115 VAC, 50/60 Hz		Description
Item Code	Model	Item Code	Model	
2170034	CCL-050B-8	2170054	CCL-050B-9	CelCulture® Incubator, 50 L, IR sensor, CO ₂ Control, Moist Heat Decon (Without Decon Pump)
2170002	CCL-170B-8	2170004	CCL-170B-9	CelCulture® Incubator, 170 L, IR sensor, CO ₂ Control, ULPA, Moist Heat Decon
2170068	CCL-170B-8-NF	2170075	CCL-170B-9-NF	CelCulture® Incubator, 170 L, IR sensor, CO ₂ Control, Moist Heat Decon (No ULPA Filter)
2170058	CCL-240B-8	2170060	CCL-240B-9	CelCulture® Incubator, 240 L, IR sensor, CO ₂ Control, ULPA, Moist Heat Decon
2170069	CCL-240B-8-NF	2170079	CCL-240B-9-NF	CelCulture® Incubator, 240 L, IR sensor, CO ₂ Control, Moist Heat Decon (No ULPA Filter)

CelCulture® CO₂ Incubators Suppressed O₂ Models with Stainless Steel Chamber

230 VAC, 50/60 Hz		115 VAC, 50/60 Hz		Description
Item Code	Model	Item Code	Model	
2170055	CCL-050T-8	2170056	CCL-050T-9	CelCulture® Incubator, 50 L, IR sensor, CO ₂ Control, Moist Heat Decon (Without Decon Pump)
2170047	CCL-170T-8	2170048	CCL-170T-9	CelCulture® Incubator, 170 L, IR sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon
2170070	CCL-170T-8-NF	2170076	CCL-170T-9-NF	CelCulture® Incubator, 170 L, IR sensor, CO ₂ & O ₂ Control, Moist Heat Decon (No ULPA Filter)
2170061	CCL-240T-8	2170062	CCL-240T-9	CelCulture® Incubator, 240 L, IR sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon
2170071	CCL-240T-8-NF	2170080	CCL-240T-9-NF	CelCulture® Incubator, 240 L, IR sensor, CO ₂ & O ₂ Control, Moist Heat Decon (No ULPA Filter)

CelCulture® CO₂ Incubator with Stainless Steel Exterior Cabinet

The Esco CelCulture® CO₂ incubator is also available with stainless steel exterior with the same superior features.

Features:

- Corrosion Resistant Surface
- Meets Pharmaceutical and Clinical Laboratory Requirements
- HPA-validated 90°C overnight moist heat decontamination cycle
- ISO Class 5 cleanliness via ULPA Filter System
- Available in 50 L, 170 L and 240 L sizes



Ordering Information

IR SENSOR MODEL WITH STAINLESS STEEL EXTERIOR CABINET		
Models	Item Code	Description
CCL-050B-8-SS	2170128	CelCulture® Incubator 50 L, IR Sensor, CO ₂ Control, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-8-SS	2170065	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-240B-8-SS	2170137	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-050B-9-SS	2170176	CelCulture® Incubator 50 L, IR sensor, CO ₂ Control, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170B-9-SS	2170177	CelCulture® Incubator 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz
CCL-240B-9-SS	2170140	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz

SUPPRESSED O ₂ MODEL WITH STAINLESS STEEL EXTERIOR CABINET		
Models	Item Code	Description
CCL-050T-8-SS	2170171	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-8-SS	2170129	CelCulture® Incubator 170 L IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-240T-8-SS	2170138	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-050T-9-SS	2170178	CelCulture® Incubator 50 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz (Without Decon Pump)
CCL-170T-9-SS	2170179	CelCulture® Incubator 170 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz
CCL-240T-9-SS	2170141	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 115 VAC, 50/60 Hz

CelCulture® CO₂ Incubator with Cooling System

Esco CelCulture® CO₂ Incubator with Integrated Cooling System provides solution for highly specialized application. The integrated cooling system allows studies of samples that requires temperature at/or below ambient temperature.

Features:

- Wider temperature range of 8°C below ambient to 60°C above ambient
- Highly efficient, environment friendly Peltier Cooling System
- Constructed using electrogalvanized with ISOCIDE™ powder-coating to eliminate 99.9% of surface bacteria within 24 hours of exposure
- Complete contamination control methods
 - 90°C Validated Moist Heat Decontamination Cycle
 - ULPA Filter
 - 0.2 µm inlet filter
- Available in 170 L and 240 L sizes



GUIDE TO MODELS, CelCulture® CO₂ Incubators with Cooling System

IR Sensor Model With Stainless Steel Chamber

230 VAC, 50/60 Hz		115 VAC, 50/60 Hz		Description
Item Code	Model	Item Code	Model	
2170101	CCL-170B-8-P	2170115	CCL-170B-9-P	CelCulture® Incubator, 170 L, IR sensor, CO ₂ control, Moist Heat Decon, Peltier System
2170116	CCL-240B-8-P	2170266	CCL-240B-9-P	CelCulture® Incubator 240 L, IR Sensor, CO ₂ Control, Moist Heat Decon, Peltier System

Suppressed O₂ Models with Stainless Steel Chamber

230 VAC, 50/60 Hz		115 VAC, 50/60 Hz		Description
Item Code	Model	Item Code	Model	
2170112	CCL-170T-8-P	2170153	CCL-170T-9-P	CelCulture® Incubator, 170 L, IR sensor, CO ₂ control, O ₂ control, Moist Heat Decon, Peltier System
2170267	CCL-240T-8-P	2170268	CCL-240T-9-P	CelCulture® Incubator 240 L, IR Sensor, CO ₂ & O ₂ Control, Moist Heat Decon, Peltier System

Highest Yield. Lowest Cost. Linearly Scalable.

A single CelCradle™ bottle replaces hundreds of petri dishes, tissue culture flasks, and dozens of spinner flasks, roller bottles, and more.



Spinner Flask
500 ml
x16



Petri Dish
180 cm²
x400



Plates
6320 cm²
x8



Flask
150 cm²
x500

1 CelCradle™ Stage accommodates
4 CelCradle™ bottles

Roller Bottles
850 cm²

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижегород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

По вопросам продаж и поддержки обращайтесь:

ADH

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

eus@nt-rt.ru || <https://escolifesciences.nt-rt.ru>